

## **Chatham Lab General IP Protocol (author: Boglarka Laczy)**

### **I. Preparation of tissue lysate (optimal cc.: 1-5 mg/mL)**

- Make lysis buffer, pipette ~ **1,5 mL** into centrifuge tubes, add ~ **100 mg** tissue powder; keep the tubes on ice;
- Homogenize the samples: 20-40 sec (clean with dd. H<sub>2</sub>O and methanol between each sample);
- Keep the homogenized samples on ice, vortex them /10-15 min;
- Centrifugation (after **45 min lysis**): **15000 x g for 15 min at 4°C**;
- Take and aliquot (350 µL) the supernatant into Eppendorf tubes;
- Store at - 80°C;

### **II. Protein A Agarose preparation**

(Protein A agarose; 10 mL packed beads containing 6 mg/mL recombinant Protein A suspended as a 50% slurry in dd H<sub>2</sub>O containing 0.01 % thiomersal; # 16-156; Upstate)

- Prior to use, **wash 2x the agarose beads** with an appropriate buffer to remove ethanol;
- Cut the tip of the pipette!
- Add PBS to Protein A slurry (1: 9); centrifuge 3 000 x g for 2 min, remove PBS;
- Add PBS to Protein A slurry (1: 10); centrifuge 3 000 x g for 2 min, remove PBS;
- Resuspend the beads in PBS and **restore to a 50 % slurry (1:1)**;

### **III. Preclearing**

- Preclear the lysate by adding **100 µL of Protein A agarose bead slurry (50%) per 1 mL of lysate (30-40 µL/sample)**;
- Incubate **at 4 °C for 60 min** with gentle agitation;
- Remove the beads by centrifugation **14 000 x g for 10 min at 4 °C**;
- **Transfer and keep the supernatant** for IP (cave beads!); beads can be washed 1-2x in PBS and collected together;

### **IV. Protein measurement**

- **Dilute the lysate** before beginning IP to roughly 1 µg/ µL total protein in a centrifuge tube with PBS;

## V. Immunoprecipitation

- **Add** 4  $\mu\text{g}$  (= 16  $\mu\text{L}$ ) of IP antibody to **500-1000  $\mu\text{g}$**  (500-1000  $\mu\text{L}$ ) the pre-cleared tissue sample lysate;

*CTR of Ab/agarose [BS]:* mock lysate (only lysis buffer without sample) plus IP Ab;

*CTR of lysate [BA]:* lysate including sample without IP Ab;

*Negative CTR [BL]:* lysate without sample and without Ab;

- Gently rock the reaction mixture **at 4°C for overnight**;
- Capture the immunocomplex by adding **100  $\mu\text{L}$**  of washed Protein A agarose bead slurry;
- Gently rock the lysate-beads mixture **at 4°C for 2hs**;
- Collect the agarose beads by pulse centrifugation: **10 000 x g for 2 min at 4°C**;
- **Drain off the supernatant**; carefully aspirate the supernatant leaving 10  $\mu\text{L}$  fluid above the pellet;
- **Wash the beads 3 x** with either ice-cold PBS (800  $\mu\text{L}$ ); Centrifuge **10 000 x g for 2 min at 4°C** and **discard the supernatant**;
- Resuspend the agarose beads in **60  $\mu\text{L}$  Laemmli sample buffer (2X)**; Do not vortex!
- Boil the beads for **5 min at 100 °C**;
- Collect the beads: by pulse centrifugation: **10 000 x g for 2 min**;
- **Take the supernatant (= IP sample)** and apply SDS-PAGE, or IP sample can be stored at -20°C (and should be reboiled prior to loading);