

UAB MRI Scanners

Before placing a patient on the scanner, ensure you have an accurate weight. No patient exceeding the limits below should be scanned. Doing so would have a high risk of causing serious scanner damage, leading to possible denial of or delay in service to patients.

Vendor	Model	Table Weight (lbs)
TKC		
GE	ARTIST	500
SIEMENS	VIDA	500
SIEMENS	SKYRA	550
GE	TWIN SPEED	350
PHILIPS	INGENIA	550
HIGHLANDS		
SIEMENS	SYMPHONY	550
SIEMENS	PRISMA	550
SIEMENS	FREEMAX	700
GE	PREMIER	550
GARDENDALE		
GE	ARTIST	500
LEEDS		
GE	ARTIST	500
ADVANCED IMAGING FACILITY (AIF)		
GE	SIGNA PET MR	500
UAB HOSPITAL NORTH PAVILION		
PHILIPS	INGENIA	550
PHILIPS	INGENIA	550
SIEMENS	SOLA	550
IMR		
SIEMENS	AREA	550