Proteomics and Thiol Modification







Aimee Landar, Ph.D.

Proteomics is the study of a protein complement in response to a stimulus





Increased protein modification in response to oxidative stress



Role of thiols in protein function and cell signaling

Detection of thiol modifications- BIAM



(Kim et al. Anal. Biochem. 283:214-221, 2000)

Assessment of Protein Modification:



liver mitochondria



Angeli's Salt





Landar and Oh, et al., FRBM 40:459, 2006

Sypro stain



veh. control



Streptavidin-HRP blot (biotin)





AS

protein stain



biotin blot



Challenges:

Quantitation of protein and thiols

Correlating thiol status with protein function

Cytochrome c as an internal standard for protein and Biotin



Cytochrome c: small (12,000 kDa), water soluble, multiple surface lysine residues.

Biotin Tagging through Lysine:



MALDI-TOF characterization of bt-cyt c

Native Cytochrome c - 12360

Matrix Adduct - 12569



Sensitivity of bt-cyt c on western blot



Internal standards for protein and Biotin



Biotinylated Cytochrome c Biotinylated Bovine serum albumin Biotinylated Myoglobin



How do we calculate the amount of biotin and protein?

Application to a 2D-Proteomic Format

Sypro Ruby stain

biotin blot





unlabeled cyt c as protein standard bt-cyt c as biotin standard



0.3µg 0.01µg



protein stain

biotin blot







Quantitation



Glutamate dehydrogenase Aldehyde dehydrogenase

mol biotin/mol protein	theoretical maximum
2.3	6
4.7	10

Relating Extent of Modification by ROS/RNS to Functional Changes



Thiol Quantitation and Proteomics

<u>Applications</u>: <u>Direct thiol measurement</u> Mitochondria Other subcellular fractions Cell lysates Tissue homogenates <u>Quantitation of different kinds of thiol</u> <u>modifications (biotin switch)</u> Reversible (-SNO, -SOH, -SSG) vs. Irreversible (SO₃H, SC_x)

The Lab

Victor Darley-Usmar JooYeun Oh Michelle Johnson Karina Ricart Anne Diers Elena Ulasova Nobuo Watanabe Jessica Gutierrez Brian Dranka Ashlee Preston Jeffrey Dubuisson

Collaborators

Stephen Barnes Amanda Isom

The Proteomic Core Facility Helen Kim, Director Landon Wilson Todd Sanderson